

DrägerService®

**Dräger**

# PMS Parts Reference Guide

Part Number: SP00217

Rev: -

Date: 31 May 2001

© 2001 Draeger Medical, Inc.

## Evita 1

[RETURN TO SERVICE PROCEDURE TABLE OF CONTENTS](#)  
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

## PMS Parts Reference Guide


This document is intended as a reference guide to assist in the replacement of routine maintenance parts.

### Evita 1

<u>Description</u>	<u>Part Number</u>	<u>Comments</u>
6 month	4114758	
	8408539	Filter 1x
	8410181	Diaphragm 2x
	8407979	Sealing Washer 1x
12 month	4114759	
	8408539	Filter 1x
	8410181	Diaphragm 2x
	8407979	Sealing Washer 1x
	8408442	Lip Seal 2x
	8407689	Lip Seal 1x
24 month	4114760	
	8408539	Filter 1x
	8410181	Diaphragm 2x
	8407979	Sealing Washer 1x
	8408442	Lip Seal 2x
	8407689	Lip Seal 1x
	2M08776	O-Ring 1x
	M20622	O-Ring 2x
	2M85714	Filter 1x
	1817582	Lithium Battery 1x
	8301856	Storage Battery 1x
	8408208	Filter 2x
	8408204	Packing 1x
	8407988	Sealing Washer 1x
72 month	<i>No Kit</i>	
	8408205	Pressure Regulator 2x


PMS Parts Reference Guide (continued)

## Evita 1




PMS Parts Reference Guide (continued)

## Expiration Valve




PMS Parts Reference Guide (continued)

## Lip Seals




PMS Parts Reference Guide (continued)

## Inspiration Block



PMS Parts Reference Guide (continued)


## Gas Supply Block






PMS Parts Reference Guide (continued)

## CPU PCB



PMS Parts Reference Guide (continued)

## Monitoring PCB



[RETURN TO SERVICE PROCEDURE TABLE OF CONTENTS](#)  
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

DrägerService®

**Dräger**

DrägerService is a division of  
Draeger Medical, Inc.  
3122 Commerce Drive  
Telford, PA 18969  
Tel: (215) 721-5402  
(800) 543-5047  
Fax: (215) 721-5784  
Web: [www.draegermedical.com](http://www.draegermedical.com)  
Printed in the U.S.A.